

	Codice prodotto - Product code - Code produit - Código producto - Erzeugniscode - Código produto - Κωδικός προϊόντος - كود المنتج
	Leggere e seguire attentamente le istruzioni per l'uso - Caution/ Read instructions carefully - Lire et suivre attentivement la notice - Leer y seguir atentamente las instrucciones de uso - Gebrauchsanweisung sorgfältig lesen und befolgen - Ler e seguir atentamente as instruções de uso - Διαβάστε και ακολουθήστε προσεκτικά τις οδηγίες χρήσης - اقرأ بدقة وحرص تعليمات الاستخدام
	Non riutilizzare - Do Not Reuse - Ne pas réutiliser - No reutilizable - Nicht wiederverwenden - Não reutilizar - Μην το επαναχρησιμοποιείτε - أداة أحادية الاستخدام، لا تستخدم من جديد
	Numero di lotto - Lot Number - Numéro de lot - Número de lote - Chargennummer - Número de lote - Αριθμός παρτίδας - رقم الدفعة (انظر العبوة / المغلف)
	Fabbricante - Manufacturer - Fabricant - Fabricante - Hersteller - Fabricante - Παραγωγός - الشركة المصنعة
	Conservare al riparo dalla luce solare - Keep away from sunlight - À conserver à l'abri de la lumière du soleil - Conservar al amparo de la luz solar - Vor Sonneneinstrahlung geschützt lagern - Guardar ao abrigo da luz solar - Κρατήστε το μακριά από ηλιακή ακτινοβολία - يحفظ بعيداً عن ضوء الشمس
	Conservare in luogo fresco ed asciutto - Keep in a cool, dry place - À conserver dans un endroit frais et sec - Conservar en un lugar fresco y seco - An einem kühlen und trockenen Ort lagern - Armazenar em local fresco e seco - Διατηρείται σε δροσερό και στεγνό περιβάλλον - يحفظ في مكان بارد وجاف
	Data di scadenza - Use By / Expiry Date - À utiliser avant le / Date de péremption - Uso por / Fecha de Caducidad - Verwendbar bis / Abluftdatum - Usar até/data de expiração - Ημερομηνία λήξεως - تاريخ انتهاء الصلاحية (انظر العبوة / المغلف)
	Prodotto conforme alla Direttiva Europea - Product complies with European Directive - Ce produit est conforme à la directive européenne - Producto conforme a la Directiva Europea - Produkt gemäß EU-Richtlinie - Produto em conformidade com a Diretiva Europeia - Προϊόν σύμφωνο με την Ευρωπαϊκή Οδηγία - الأوروبي الاتحادي توافق مع يتوافق منتج
	Leggere attentamente le istruzioni per l'uso - Please read instructions carefully - Lire attentivement la notice - Leer atentamente las instrucciones de uso - Gebrauchsanweisung sorgfältig lesen - Ler atentamente as instruções de uso - Διαβάστε προσεκτικά τις οδηγίες χρήσης - اقرأ بدقة وحرص تعليمات الاستخدام
	Data di fabbricazione - Date of Manufacturing - Date de fabrication - Fecha de fabricación - Herstellungsdatum - Data de fabrico - Ημερομηνία παραγωγής - تاريخ التصنيع



PROFESSIONAL MEDICAL PRODUCTS

Gima S.p.A.
Via Marconi, 1 - 20060 Gessate (MI) Italy
gima@gimaitaly.com - export@gimaitaly.com
www.gimaitaly.com

**MASCHERINA SALVAVITA CPR
CPR LIFE MASK
MASQUE DE POCHE (RCP)
NOTFALLBEATMUNGSMASKE
MÁSCARA SALVAVIDA
MÁSCARA SALVA-VIDA
ΜΑΣΚΑ ΔΙΑΣΩΣΗΣ**

قناع خاص لعملية الإنقاذ

Manuale d'uso e manutenzione
Use and maintenance book
Instructions de fonctionnement
et entretien
Betriebs- und wartungs-
anweisungen
Manual de uso y mantenimiento
Manual de uso e manutenção
Εγχειρίδιο χρήσης
και συντήρησης
دليل الإستعمال والرعاية



ATTENZIONE: Gli operatori devono leggere e capire completamente questo manuale prima di utilizzare il prodotto.

ATTENTION: The operators must carefully read and completely understand the present manual before using the product.

AVIS: Les opérateurs doivent lire et bien comprendre ce manuel avant d'utiliser le produit.

ACHTUNG: Diese Anleitung muss vor dem Einsatz des Produkts aufmerksam gelesen und vollständig verstanden werden.

ATENCIÓN: Los operadores tienen que leer y entender completamente este manual antes de utilizar el producto.

ATENÇÃO: Os operadores devem ler e entender completamente este manual antes de usar o produto.

ΠΡΟΣΟΧΗ: Οι χειριστές αυτού του προϊόντος πρέπει να διαβάσουν και να καταλάβουν πλήρως τις οδηγίες του εγχειριδίου πριν από την χρήση του.

الحنن: على العمال قراءة وفهم هذا الدليل بكامله قبل البدء باستخدام المنتج.

M34218-M-Rev3-01.19


REF 34218

Gima S.p.A.
Via Marconi, 1
20060 Gessate (MI) Italy
Made in China







A




B



C



D



CPR LIFE MASK

Mouth-to-mouth respiration adult/children

This product was designed only for medical or paramedical staff specialized in resuscitation techniques.

- (A) Lay patient on back with a pillow or similar object under the shoulders to keep the patient straight and free the mouth.
- (B) Place mask with the tighter part over the nose and the anti-reflux valve (OUT) towards the chest of the patient and secure it to the head using the elastic.
- (C) If oxygen should be necessary connect the tube to the specific connection. In this case verify that there are no flames, oils or grease nearby.
- (D) The mask is a disposable device, it cannot be put in an autoclave. After use, throw away and dispose of according to the laws in force.



CAUTION: *Gima S.p.A. cannot be held responsible for malfunctioning or damages resulting from tampering with the device, abuse, incorrect storage or bad environmental conditions, or for incorrect use on behalf of unqualified or unauthorized personnel.*